

AXEL

US Biological (antibodymab)

85-8564-22 Phosphoglycerate Mutase 1 (PGAM1, PGAM-1, PGAMA, BPG-dependent PGAM1, BPG-dependent PGAM-1, OTTHUMP00000059414, PGAM-B, Phosphoglycerate Mutase A, Phosphoglycerate Mutase Isozyme B) (MaxLight 650) 100ul 131287-ML650

USBiological
Life Sciences

アズワン品番	85-8564-22
型番	131287-ML650
入り数	1個
標準価格	-円

USBiological
Life Sciences

仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

650.

Specificity: Recognizes human PGAM1.

Clone # USB:2G1-A6

Isotype: IgG1,k

Calc Applications Abbrev: FLISA IP WB

Calc Crossreactivity: Hu

Immunogen: Full length recombinant corresponding to aa1-255 from human PGAM1 (AAH11678) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1year

EU Commodity Code: 30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp