

US Biological (antibodymab)

85-8564-24 Phospholipase C, gamma 1 (Phospholipase C-gamma-1, PLCgamma1, PLC-gamma-1, PLCg1, NCKAP3, 1-Phosphatidylinositol-4,5-bisphosphate Phosphodiesterase gamma-1, Phospholipase C-II, PLC-II, Phosphoinositide Phospholipase C-gamma-1, PLC1, PLC148, PLC-148) (MaxLight 650) 100ul 131290-ML650



アズワン品番	85-8564-24
型番	131290-ML650
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

650.

Specificity: Recognizes human PLCG1.

Clone # USB: 2A2

Isotype: IgG1, k

Calc Applications Abbrev: FLISA IF WB

Calc Crossreactivity: Hu

Immunogen: Partial recombinant corresponding to aa1192-1292 from human PLCG1 (NP_002651) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1 year

EU Commodity Code: 30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp