

US Biological (antibodymab)

85-8564-25 PHOX2A (Paired Mesoderm Homeobox Protein 2A, ARIX1 Homeodomain Protein, Aristaless Homeobox Protein Homolog, Paired-like Homeobox 2A, ARIX, PMX2A) (MaxLight 650) 100ul 131292-ML650



アズワン品番	85-8564-25
型番	131292-ML650
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

tTM650.

Specificity: Recognizes human PHOX2A. Species Crossreactivity: rat.

Clone # USB:4F6

Isotype:IgG2a,k

Calc Applications Abbrev:FLISA IF WB

Calc Crossreactivity:Hu Rt

Immunogen:Partial recombinant corresponding to aa1-90 from human PHOX2A (NP_005160) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life:1year

EU Commodity Code:30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp