

US Biological (antibodymab)

85-8564-28 PI15 (Peptidase Inhibitor 15, PI-15, 25k D Trypsin Inhibitor, p25TI, Cysteine-rich Secretory Protein 8, CRISP-8, SugarCrisp, CRISP8, P25TI, DKFZp686F0366, P24TI) (MaxLight 650) 100ul 131298-ML650



アズワン品番	85-8564-28
型番	131298-ML650
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

650.

Specificity: Recognizes human PI15.

Clone # USB:3B5

Isotype:IgG1,k

Calc Applications Abbrev:FLISA IHC WB

Calc Crossreactivity:Hu

Immunogen:Partial recombinant corresponding to aa18-118 from human PI15 (NP_056970) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life:1year

EU Commodity Code:30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp