

US Biological (antibodymab)

85-8564-29 PI4KB (PIK4CB, Phosphatidylinositol 4-kinase beta, PI4K-beta, PI4Kbeta, PtdIns 4-kinase beta, NPIK, PI4K92) (MaxLight 650) 100ul 131304-ML650



アズワン品番	85-8564-29
型番	131304-ML650
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight 650.

Specificity: Recognizes human PIK4CB.

Clone # USB: 3B1

Isotype: IgG2b,k

Calc Applications Abbrev: FLISA WB

Calc Crossreactivity: Hu

Immunogen: Partial recombinant corresponding to aa731-829 from human PIK4CB (NP_002642) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1year

EU Commodity Code: 30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp