



US Biological (antibodymab)

85-8564-35 PICK1 (Protein Interacting with PRKCA 1, Protein Interacting with C Kinase 1, Protein that Interacts with C Kinase 1, PICK-1, MGC15204, OTTH UMP00000028509, Protein Kinase C alpha Binding Protein, PRKCA Binding Protein, Prkcabp) (MaxLight 650) 100ul 131319-ML650



アズワン品番	85-8564-35
型番	131319-ML650
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

650.

Specificity: Recognizes human PRKCABP. Species Crossreactivity: rat.

Clone # USB: 3G5

Isotype: IgG2a,k

Calc Applications Abbrev: FLISA WB

Calc Crossreactivity: Hu Rt

Immunogen: Partial recombinant corresponding to aa266-375 from PRKCABP (NP_036539) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1 year

EU Commodity Code: 30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp