

US Biological (antibodymab)

85-8564-40 PIGO (GPI Ethanolamine Phosphate Transferase 3, Phosphatidylinositol-glycan Biosynthesis Class O Protein, PIG-O, UNQ632/PRO1249, DKFZp434M222, FLJ00135, MGC20536, MGC3079, RP11-182N22.4) (MaxLight 650) 100ul 131326-ML650



アズワン品番	85-8564-40
型番	131326-ML650
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

650.

Specificity: Recognizes human PIGO.

Clone # USB:2A4-1C3

Isotype:IgG2b,k

Calc Applications Abbrev:FLISA WB

Calc Crossreactivity:Hu

Immunogen:Full length recombinant corresponding to aa1-455 from human PIGO (AH01030) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life:1year

EU Commodity Code:30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp