

US Biological (antibodymab)

85-8564-41 PIGP (Phosphatidylinositol N-acetylglucosaminyltransferase Subunit P, Down Syndrome Critical Region Protein 5, Down Syndrome Critical Region Protein C, Phosphatidylinositol-glycan Biosynthesis Class P Protein, PIG-P, DCRC, DSCR5, DSCRC, NPD010) (MaxLight 650) 100ul 131327-ML650



アズワン品番	85-8564-41
型番	131327-ML650
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

650.

Specificity: Recognizes human PIGP.

Clone # USB: 2E7

Isotype: IgG2a,k

Calc Applications Abbrev: FLISA WB

Calc Crossreactivity: Hu

Immunogen: Partial recombinant corresponding to aa77-135 from human PIGP (NP_710149) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1 year

EU Commodity Code: 30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp