

US Biological (antibodymab)

85-8564-48 PIK3C2B (Phosphatidylinositol-4-phosphate 3-kinase C2 Domain-containing Subunit beta, PI3K-C2-beta, Phosphoinositide 3-kinase-C2-beta, PIK3C2B, PtdIns-3-kinase C2 Subunit beta, C2-PI3K, DKFZp686G16234) (MaxLight 650) 100ul 131338-ML650



アズワン品番	85-8564-48
型番	131338-ML650
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

650.

Specificity: Recognizes human PIK3C2B.

Clone # USB: 3E5

Isotype: IgG2a,k

Calc Applications Abbrev: FLISA IF WB

Calc Crossreactivity: Hu

Immunogen: Partial recombinant corresponding to aa1-111 from human PIK3C2B (NP_002637) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1 year

EU Commodity Code: 30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp