

US Biological (antibodymab)

85-8564-53 PIK3R1 (GRB1, Phosphatidylinositol 3-kinase Regulatory Subunit alpha, Phosphatidylinositol 3-kinase 85kD Regulatory Subunit alpha) (MaxLight 650) 100ul 131345-ML650



アズワン品番	85-8564-53
型番	131345-ML650
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

650.

Specificity: Recognizes human PIK3R1. Species Crossreactivity: mouse and rat.

Clone # USB: 3A10

Isotype: IgG1, k

Calc Applications Abbrev: FLISA IP WB

Calc Crossreactivity: Hu Mo Rt

Immunogen: Full length recombinant corresponding to aa1-454 from human PIK3R1 (AH30815) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1year

EU Commodity Code: 30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp