

US Biological (antibodymab)

85-8564-59 PIN1 (Peptidyl-prolyl Cis-trans Isomerase NIMA-interacting 1, Peptidyl-prolyl Cis-trans Isomerase Pin1, PPlase Pin1, Rotamase Pin1) (MaxLight 650) 100ul 131353-ML650



アズワン品番	85-8564-59
型番	131353-ML650
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight 650.

Specificity: Recognizes human PIN1.

Clone # USB: 2F2

Isotype: IgG2b,k

Calc Applications Abbrev: FLISA WB

Calc Crossreactivity: Hu

Immunogen: Partial recombinant corresponding to aa64-163 from human PIN1 (AAH02899) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1year

EU Commodity Code: 30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp