

US Biological (antibodymab)

85-8564-60 PIP4K2A (Phosphatidylinositol-5-phosphate 4-kinase Type-2 alpha, 1-phosphatidylinositol-5-phosphate 4-kinase 2-alpha, Diphosphoinositide Kinase 2-alpha, PIP5KIII, Phosphatidylinositol-5-phosphate 4-kinase Type II alpha, PI(5)P 4-kinase Type II alpha, PIP4KII-alpha, PtdIns(4)P-5-kinase B Isoform, PtdIns(4)P-5-kinase C Isoform, PtdIns(5)P-4-kinase Isoform 2-alpha, PIP5K2, PIP5K2A) (MaxLight 650) ...
131357-ML650

USBiological
Life Sciences

アズワン品番	85-8564-60
型番	131357-ML650
入り数	1個
標準価格	-円

USBiological
Life Sciences

仕様

商品名: PIP4K2A (Phosphatidylinositol-5-phosphate 4-kinase Type-2 alpha, 1-phosphatidylinositol-5-phosphate 4-kinase 2-alpha, Diphosphoinositide Kinase 2-alpha, PIP5KIII,

Phosphatidylinositol-5-phosphate 4-kinase Type II alpha, PI(5)P 4-kinase Type II alpha, PIP4KI-alpha, PtdIns(4)P-5-kinase B Isoform, PtdIns(4)P-5-kinase C Isoform, PtdIns(5)P-4-kinase Isoform 2-alpha, PIP5K2, PIP5K2A) (MaxLight 650) 100ul

Size:100ul

Host:mouse

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by Protein A affinity chromatography.

Form:Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight™650.

Specificity:Recognizes human PIP5K2A.

Clone # USB:3A3

Isotype:IgG1,k

Calc Applications Abbrev:FLISA IHC WB

Calc Crossreactivity:Hu

Immunogen:Partial recombinant corresponding to aa304-366 from human PIP5K2A (NP_005019) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life:1year

EU Commodity Code:30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp