

US Biological (antibodymab)

85-8564-61 PIP5K1C (KIAA0589, Phosphatidylinositol 4-phosphate 5-kinase Type-1 gamma, Phosphatidylinositol 4-phosphate 5-kinase Type I gamma) (MaxLight 650) 100ul 131363-ML650



|        |              |
|--------|--------------|
| アズワン品番 | 85-8564-61   |
| 型番     | 131363-ML650 |
| 入り数    | 1個           |
| 標準価格   | -円           |



## 仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight 650.

Specificity: Recognizes human PIP5K1C.

Clone # USB: 7D11

Isotype: IgG2b,k

Calc Applications Abbrev: FLISA WB

Calc Crossreactivity: Hu

Immunogen: Partial recombinant corresponding to aa561-668 from PIP5K1C (NP\_036530) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1year

EU Commodity Code: 30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp