

US Biological (antibodymab)

85-8564-62 PIP5K3 (PIKFYVE, 1-phosphatidylinositol 3-phosphate 5-kinase, Phosphatidylinositol 3-phosphate 5-kinase, FYVE Finger-containing Phosphoinositide Kinase, PIKfyve, Phosphatidylinositol 3-phosphate 5-kinase Type III, PIPkin-III, Type III PIP Kinase, KIAA0981, FLJ37746, MGC40423) (MaxLight 650) 100ul 131366-ML650

USBiological
Life Sciences

アズワン品番	85-8564-62
型番	131366-ML650
入り数	1個
標準価格	-円

USBiological
Life Sciences

仕様

Size:100ul

Host:mouse

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

650.

Specificity: Recognizes human PIP5K3.

Clone # USB: 6C7

Isotype: IgG2a,k

Calc Applications Abbrev: FLISA IF IHC WB

Calc Crossreactivity: Hu

Immunogen: Partial recombinant protein corresponding to aa342-452 from human PIP5K3 (NP_689884) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1 year

EU Commodity Code: 30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp