

US Biological (antibodymab)

85-8564-69 PITX1 (Pituitary Homeobox 1, Paired-like Homeodomain Transcription Factor 1, Homeobox Protein PITX1, Hindlimb-expressed Homeobox Protein Backfoot, PTX1, BFT) (MaxLight 650) 100ul 131377-ML650



アズワン品番	85-8564-69
型番	131377-ML650
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

650.

Specificity: Recognizes human PITX1.

Clone # USB:5G4

Isotype:IgG2a,k

Calc Applications Abbrev:FLISA IP WB

Calc Crossreactivity:Hu

Immunogen:Partial recombinant corresponding to aa225-314 from human PITX1 (NP_002644) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life:1year

EU Commodity Code:30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp