

US Biological (antibodymab)

85-8564-71 PJA1 (E3 Ubiquitin-protein Ligase Praja-1, PRAJA1, Ring Finger Protein 70, RNF70) (MaxLight 650) 100ul 131383-ML650

USBiological
Life Sciences

アズワン品番	85-8564-71
型番	131383-ML650
入り数	1個
標準価格	-円

USBiological
Life Sciences

仕様

- Size: 100ul
- Host: mouse
- Source Antibody: human
- Grade: Affinity Purified
- Purity: Purified by Protein A affinity chromatography.
- Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight™650.
- Specificity: Recognizes human PJA1.
- Clone # USB: 3E10
- Isotype: IgG2a,k
- Calc Applications Abbrev: FLISA WB
- Calc Crossreactivity: Hu
- Immunogen: Partial recombinant corresponding to aa81-181 from human PJA1 (NP_660095) with GST tag. MW of the GST tag alone is 26kD.
- Shelf Life: 1year
- EU Commodity Code: 30021010
- この商品は法規制を確認していません。(法規制によって販売できない場合もございます)
- 製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。
- 【試薬に関するお問合せ】
- アズワン株式会社 試薬・プロセス材料グループ
- TEL: 06-6447-8641
- FAX: 06-6447-8642
- E-mail: Reagents@so.as-1.co.jp