

US Biological (antibodymab)

85-8564-73 PKD2 (Polycystin-2, Autosomal Dominant Polycystic Kidney Disease Type II Protein, Polycystic Kidney Disease 2 Protein, Polycystin, R48321) (MaxLight 650) 100ul 131388-ML650



アズワン品番	85-8564-73
型番	131388-ML650
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight 650.

Specificity: Recognizes human PKD2.

Clone # USB: 4C8

Isotype: IgG2a,k

Calc Applications Abbrev: FLISA IF WB

Calc Crossreactivity: Hu

Immunogen: Partial recombinant corresponding to aa261-361 from human PKD2 (NP_000288) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1year

Swiss Prot Number: Q13563

EU Commodity Code: 30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp