

US Biological (antibodymab)

85-8564-78 PKN1 (Serine/Threonine-protein Kinase N1, Protease-activated Kinase 1, PAK-1, Protein Kinase C-like 1, Protein Kinase C-like PKN, Protein Kinase PKN-alpha, Protein-kinase C-related Kinase 1, Serine-threonine Protein Kinase N, PKN1, PAK1, PRK1, PRKCL1, MGC46204) (MaxLight 650) 100ul 131401-ML650

**USBiological**  
Life Sciences

|        |              |
|--------|--------------|
| アズワン品番 | 85-8564-78   |
| 型番     | 131401-ML650 |
| 入り数    | 1個           |
| 標準価格   | -円           |

**USBiological**  
Life Sciences

## 仕様

Size:100ul

Host:mouse

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

650.

Specificity: Recognizes human PKN1.

Clone # USB: 1B10

Isotype: IgG2a,k

Calc Applications Abbrev: FLISA IF WB

Calc Crossreactivity: Hu

Immunogen: Partial recombinant corresponding to aa462-615, from human PKN1 (AA H40061) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1year

EU Commodity Code: 30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp