

US Biological (antibodymab)

85-8564-96 PLEKHC1 (FERMT2, KIND2, MIG2, Fermitin Family Homolog 2, Kindlin-2, Mitogen-inducible Gene 2 Protein, Pleckstrin Homology Domain-containing Family C Member 1, DKFZp686G11125, FLJ34213, FLJ44462, UNC112, UNC112B, Mig-2) (MaxLight 650) 100ul 131455-ML650

USBiological
Life Sciences

アズワン品番	85-8564-96
型番	131455-ML650
入り数	1個
標準価格	-円

USBiological
Life Sciences

仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

650.

Specificity: Recognizes human PLEKHC1.

Clone # USB:2G11

Isotype: IgG2a,k

Calc Applications Abbrev: FLISA IF WB

Calc Crossreactivity: Hu

Immunogen: Full length recombinant corresponding to aa1-681 from human PLEKHC1 (AAH17327) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1year

EU Commodity Code: 30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp