

US Biological (antibodymab)

85-8564-99 PLEKHM1 (Pleckstrin Homology Domain-containing Family M Member 1, PH Domain-containing Family M Member 1, 162kD Adapter Protein, A P162, KIAA0356, OPTB6) (MaxLight 650) 100ul 131460-ML650



アズワン品番	85-8564-99
型番	131460-ML650
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

650.

Specificity: Recognizes human PLEKHM1.

Clone # USB:1C9

Isotype:IgG2b,k

Calc Applications Abbrev:FLISA

Calc Crossreactivity:Hu

Immunogen:Partial recombinant corresponding to aa957-1056 from PLEKHM1 (AAH 64361) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life:1year

EU Commodity Code:30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp