



G-Biosciences

RED 660 Protein Assay, 500 Assays/ 2,500 Micro-assays

特徴

Protein estimation of samples in aqueous buffers

Estimation in presence of detergents



シリーズ共通仕様

RAP ID	Single reagent assays	VERSATILE	Compatible with higher range of detergents and reducing agents
LINEAR	Perfect linear standard curves compared to other protein assays	容量	500 Assays/ 2,500 Micro-assays

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
--------	--------	-----	----	----	------------------------------	----



89-5259-68

RED 660 Protein Assay, 500 Assays
/ 2,500 Micro-assays

786-676

1個

63,000円



89-5259-69

RED 660 Protein Assay with Non
Animal Protein Standard, 500 Assays/
2,500 Micro-assays

786-899

1個

63,000円