



G-Biosciences

RED 660 Protein Assay with Non Animal Protein Standard, 500 Assays/ 2,500 Micro-assays

特徴

- Protein estimation of samples in aqueous buffers
- Estimation in presence of detergents



シリーズ共通仕様

RAP ID	Single reagent assays	VERSATILE	Compatible with higher range of detergents and reducing agents
LINEAR	Perfect linear standard curves compared to other protein assays	容量	500 Assays/ 2,500 Micro-assays

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
	89-5259-68	RED 660 Protein Assay, 500 Assays / 2,500 Micro-assays	786-676	1個	63,000円	
	89-5259-69	RED 660 Protein Assay with Non Animal Protein Standard, 500 Assays / 2,500 Micro-assays	786-899	1個	63,000円	