

G-Biosciences

RIPA Lysis & Extraction Buffer, 100mL

特徴

A complete lysis buffer for the release of cytoplasmic, membrane and nuclear proteins from adherent and suspension cultured mammalian cells. The RIPA lysis buffer is fully compatible with many applications, including reporter assays, protein assays, immunoassays and other protein purification techniques.



商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
	89-5259-77	RIPA Lysis & Extraction Buffer, 100mL	786-489	1個	27,000円	
	89-5259-78	RIPA Lysis & Extraction Buffer, 500mL	786-490	1個	59,000円	
	89-5259-79	RIPA Lysis & Extraction Buffer, 1L	786-723	1個	97,000円	
	89-5259-80	RIPA Lysis & Extraction Buffer, 1 Gallon	786-746	1個	297,000円	