



G-Biosciences

RUBEO Fluorescent Protein Stain, 1L

特徴

For staining of protein electrophoresis gels

Mass spectrometry analysis



シリーズ共通仕様

Fluorescent	Red emission at 605nm.	High Sensitivity	Easily detecting levels of protein.
Minimal protein protein variation	Allows quantitative comparison of proteins.	Compatible	Ideal for subsequent analysis, including mass spectrometry, microsequencing and immunostaining.
Versatile	Stain native, denaturing, 1D, 2D and 3D gels.	Suitability	Visualize with the simplest UV transilluminators to the complex laser scanners

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
--------	--------	-----	----	----	------------------------------	----



89-5259-85

RUBEO Fluorescent Protein Stain,
200mL

786-644

1個

22,000円

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
	89-5259-86	RUBEO Fluorescent Protein Stain, 1L	786-645	1個	73,000円	