

シグマアルドリッチジャパン(Sigma-Aldrich)

89-8564-23 Trypsin Inhibitor from Glycine max (soybean), Suitable for manufacturing of diagnostic kits and reagents, suitable for cell culture SRE0044-250MG CAS No:9035-81-8 SRE0044

[Sigma Aldrich Japan(Sigma-Aldrich)] Trypsin Inhibitor from Glycine max (soybean), Suitable for manufacturing of diagnostic kits and reagents, suitable for cell culture SRE0044-250MG CAS No:9035-81-8

アズワン品番	89-8564-23
型番	SRE0044
入り数	1個
標準価格	-円

仕様

別名:Trypsin inhibitor from Glycine max (soybean); SBTI

サイズ:250MG

CAS No:9035-81-8

バルクサイズご要望の場合はお問い合わせください。

英語SDSの検索はこちら。[英語SDS] <http://www.sigmaaldrich.com/safety-center.html>

※国内外の在庫状況により納期が表記と異なる事がございますことご了承ください。

CAS No :9035-81-8